

SIMULTANEOUS SUCTION AND IRRIGATION FOR EASY, COMFORTABLE SALINE IRRIGATION

Performing saline irrigation for your patients without causing them discomfort or creating a mess in the Ambulatory Centre is now a reality. The Cyclone™ sinonasal suction and irrigation system provides **SIMULTANEOUS** suction and saline irrigation in one **SINGLE-HANDED** instrument. It's unique dual-action system is designed to **SUCTION MORE THAN IT IRRIGATES**, minimising pooling and drainage of effluent to improve your patient's comfort.

FUNCTIONALITY THAT'S AHEAD OF THE CURVE

MALLEABLE TIP
ALLOWS DEVICE TO
BE MODIFIED FOR
USE IN ALL SINUSES
AS APPROPRIATE

ADJUSTABLE SUCTION
TUBE PROVIDES
VARIABLE POSITIONING
TO CAPTURE ALL
EFFLUENT REGARDLESS
OF THE ANATOMY

6 PORT
IRRIGATION TIP
PROVIDES A THOROUGH
RINSE IN ALL DIRECTIONS

PEELABLE EXTENSION
TUBES ALLOW
CUSTOMISATION
IN SET-UP



THE UNIQUE SINONASAL DUAL-ACTION SUCTION AND SALINE IRRIGATION SYSTEM



1 DEVICE

A single device provides simultaneous suction and saline irrigation. Reshape the tip to fit all of an individual patient's appropriate sinuses.

1 HAND

Single-handed placement makes the Cyclone system easy to use.

1 PACKAGE

All of the components you need come in a single package.

Stop avoiding saline irrigation in the Ambulatory Centre. With the easy-to-use Cyclone sinonasal suction and irrigation system you can comfortably and effectively perform saline irrigation in your day cases. Contact your Entellus Medical Representative or visit EntellusMedical.com/UK/Cyclone.

To order, please call **+1 800 3586046**

CATALOG NO.	DESCRIPTION
CC-100-I	Cyclone sinonasal suction and irrigation system



Entellus Medical, Inc.
3600 Holly Lane North, Suite 40
Plymouth, MN 55447

☎ 866.620.7615
☎ 763.463.1599
www.EntellusMedical.com
www.SinusSurgeryOptions.com

Products referenced are approved for CE Mark. Products may not be available in all countries.
CAUTION: Federal (USA) law restricts these devices to sale by or on the order of a physician.
Please see Instructions for Use (IFU) for a complete listing of warnings, precautions, and adverse events.
Indication for Use: To provide a means to irrigate and suction within sinonasal spaces.
Entellus and Cyclone are trademarks of Entellus Medical, Inc.

3114-036 rA 06/2016

© 2016 Entellus Medical, Inc.